



LABORATORY
EQUIPMENT

stirrers

WWW.FALCINSTRUMENTS.IT

ELECTROMAGNETIC STIRRERS

Small volume magnetic stirrers without heating



F20

-  Useful in the microbiological and biochemical fields
-  The stainless steel plate facilitates cleaning and sterilization
-  Electronic speed control DC motor

Technical data		F20
Stainless steel working surface LxW	mm	110x110
Max stirring capacity with H ₂ O	lt	2
Variable speed	Rpm	100-1100
Hole for rod	M 8	No
Motor power	W	10
Dimensions LxWxH	mm	110x115x60
Weight	kg	0,7
Stirring bar suggested Ø	mm	6x20
Code		601.0101.20

V/HZ 220/240-50/60

SAFETY CLASS 0

PROTECTION CLASS IP 41

STIRRER

Stirrer without heating - F20

-  Useful in the microbiological and biochemical fields
-  The stainless steel plate facilitates cleaning and sterilization
-  Electronic speed control DC motor

Electromagnetic stirrer F20 is a model without heating, suitable for continuous applications. It is a friendly product suggested for microbiological and biochemical fields.

Variable speed : 100 - 1100 Rpm
Max stirring capacity : 2 Lt.



[illegible]



FALC INSTRUMENTS s.r.l.

Via G. M. Compagnoni, 2
24047 Treviglio (BG) - Italy
+39 0363 304660 | falc@falcinstruments.it

WWW.FALCINSTRUMENTS.IT