



LABORATORY
EQUIPMENT

stirrers

WWW.FALCINSTRUMENTS.IT

ELECTROMAGNETIC STIRRERS

Multi-position stirrers



Digital timer



Plate in stainless steel (304). The stainless steel plate facilitates cleaning and sterilization.



Case in stainless steel coated, painted in oven with epoxy powder antacid



Displays simultaneous three function: speed, time and countdown



Digital speed control
DC motor



F4 - F8
WITHOUT HEATING

Technical data

		F4	F8
Working plate	mm	205x205	205x410
Distance from stirring points	mm	100	100
Max stirring capacity with H ₂ O	lt	1x4	1x8
Timer	hh:mm	∞ o 99:59	∞ o 99:59
Variable speed	rpm	80÷1200	80÷1200
Motor Power	W	25	25
External dimensions LxWxH	mm	205x280x70	205x485x70
Weight	kg	2,0	4,0
Stirring bar suggested Ø	mm	6x20	6x20
Code		601.0101.44	601.0101.49

V/HZ 220/240-50/60

SAFETY CLASS 0

PROTECTION CLASS IP 40

[illegible]



FALC INSTRUMENTS s.r.l.

Via G. M. Compagnoni, 2
24047 Treviglio (BG) - Italy
+39 0363 304660 | falc@falcinstruments.it

WWW.FALCINSTRUMENTS.IT