



LABORATORY
EQUIPMENT

stirrers

F50K - F90K





WWW.FALCINSTRUMENTS.IT

ELECTROMAGNETIC STIRRERS

Small volume magnetic stirrers without heating

High volume magnetic stirrers without heating



-  Ideal during pre-production process and on a vast scale
-  Stainless steel plate and metal structure painted with anti-acid epoxy powder coating
-  Digital speed control DC motor
-  Digital timer

Technical specifications		F50K	F90K
Stainless steel working surface LxW	mm	300x300	400x400
Max stirring capacity with H ₂ O	lt	50	90
Variable speed	Rpm	150-1500	150-1200
Timer		Up to 9 hours	Up to 9 hours
Max stirring bar lenght	mm	110	155
Stirring bar suggested		601.0125.71	601.0125.72
Motor power	W	15	20
Dimensions LxWxH	mm	300x410x70	400x510x110
Weight	kg	6	10
Code		601.0101.55	601.0101.56

V/HZ 110-230/50-60

SAFETY CLASS 0

PROTECTION CLASS IP 42

2022

● lincolmostudio.it



FALC INSTRUMENTS s.r.l.

Via G. M. Compagnoni, 2
24047 Treviglio (BG) - Italy
+39 0363 304660 | falc@falcinstruments.it

WWW.FALCINSTRUMENTS.IT