

ANALYTICAL BALANCES

PISA MODELS

PISA BALANCES, for their advanced performance and technical characteristics are now the more competitive on the market , for their high reliability and design they establish the best opportunities in the laboratory balance field. Small, yet compact, and with high quality technology , PISA balances feature a wide range of capacities and resolve all weighing problems in the laboratory. Possibility of automatic calibration with external weight, multifunctional green keypad (grams, carats, ounces, piece-counter) subtractive tare memory for all models, whilst T models are not provided with multifunctional keypad. FA Series are complete of weight, tweezers and brush. All balances are supplied with instruction manual.


SERIES FA - A



Series		FA 604 A	FA 1004 A	FA 1104 A	FA 1604 A	FA 2004 A	FA 2104 A
Capacity	g	0 - 60	0 - 100	0 - 110	0 - 160	0 - 200	0 - 210
Precision	g	0,0001	0,0001	0,0001	0,0001	0,0001	0,0001
Plate Ø	mm	80	80	80	80	80	80
Output data	Rs232C	yes	yes	yes	yes	yes	yes
Dimensions LxWxH	mm	200x320x290	200x320x290	200x320x290	200x320x290	200x320x290	200x320x290
Power supply	V/Hz	220/50	220/50	220/50	220/50	220/50	220/50
Item		630.0060.04	630.0100.04	630.0110.04	630.0160.04	630.0200.04	630.0210.04
Calibration weight	g	60	100	100	100	200	200


SERIES FA - AS DUAL RANGE

Series		FA 1604 AS	FA 2104 AS
Capacity	g	0 -30 / 160	0 -60 / 210
Precision	g	0,0001 / 0,001	0,0001 / 0,001
Plate Ø	mm	80	80
Output data	Rs232C	yes	yes
Dimensions LxWxH	mm	200x320x290	200x320x290
Power supply	V/Hz	220/50	220/50
Item		630.1160.04	630.1200.04
Calibration weight	g	100	200




ANALYTICAL & TECHNICAL BALANCES

JA - A SERIES



Series		JA 1003 A	JA 1203 A	JA 2003 A	JA 3003 A	JA 4103 A	JA 5003 A
Capacity	g	0 - 100	0 - 120	0 - 200	0 - 300	0 - 410	0 - 500
Precision	g	0,001	0,001	0,001	0,001	0,001	0,001
Plate Ø	mm	80	80	80	110	110	110
Output data	Rs232C	yes	yes	yes	yes	yes	yes
Dimensions LxWxH	mm	200x320x290	200x320x290	200x320x290	200x320x290	200x320x290	200x320x290
Power supply	V/Hz	220/50	220/50	220/50	220/50	220/50	220/50
Item		630.1200.03	630.1201.03	630.1202.03	630.1300.03	630.1400.03	630.1500.03
Calibration weight	g	100	100	200	200	200	200

TECHNICAL BALANCES JA SERIES



Series		JA 5002	JA 6102	JA 8002	JA 10002
Capacity	g	0 - 500	0 - 610	0 - 800	0 - 1000
Precision	g	0,01	0,01	0,01	0,01
Plate Ø	mm	160	160	160	160
Output data	Rs232C	yes	yes	yes	yes
Dimensions LxWxH	mm	200x320x80	200x320x80	200x320x80	200x320x80
Power supply	V/Hz	220/50	220/50	220/50	220/50
Item		630.2500.02	630.2600.02	630.2800.02	630.3100.02
Calibration weight	g	200	200	500	1000

Series		JA 12002	JA 21002	JA 31002
Capacity	g	0 - 1200	0 - 2100	0-3100
Precision	g	0,01	0,01	0,01
Plate Ø	mm	160	160	160
Output data	Rs232C	yes	yes	yes
Dimensioni LxWxH	mm	200x320x80	200x320x80	200x320x80
Power supply	V/Hz	220/50	220/50	220/50
Item		630.3101.02	630.3102.02	630.3103.02
Calibration weight	g	1000	1000	2000

TECHNICAL BALANCES

T SERIES



Series		T 200	T 500	T 1000	T 2000	T 3000	T 5000
Capacity	g	0 - 200	0 - 500	0 - 1000	0 - 2000	0 - 3000	0 - 5000
Precision	g	0,1	0,1	0,1	0,5	0,5	1,0
Plate	mm	Ø 140	Ø 140	Ø 140	LxP 190x160	LxP 190x160	LxP 190x160
Output data	Rs232C	no	no	no	no	no	no
Dimensions LxWxH	mm	195x215x62	195x215x62	195x215x62	195x215x62	195x215x62	195x215x62
Power supply	V/Hz	220/50	220/50	220/50	220/50	220/50	220/50
Item		630.3020.10	630.3050.10	630.3100.10	630.3200.50	630.3300.50	630.3501.00
Calibration weight	g	200	500	1000	2000	2000	5000

PORTABLE BALANCES

PT SERIES

Series		PT 102	PT 202	PT 302
Capacity	g	0 - 100	0 - 200	0 - 300
Precision	g	0,01	0,01	0,01
Plate Ø	mm	110	110	110
Output data	Rs232C	yes	yes	yes
Dimensions LxWxH	mm	175x220x70	175x220x70	175x220x70
Power supply		Battery DC9V	Battery DC9V	Battery DC9V
Item		630.4100.02	630.4200.02	630.4300.02
Calibration weight	g	100	200	200



Series		PT 401	PT 601	PT 801	PT 1001	PT 1201	PT 2101	PT 3101
Capacity	g	0 - 400	0 - 600	0 - 800	0 - 1000	0 - 1200	0 - 2100	0 - 3100
Precision	g	0,1	0,1	0,1	0,1	0,1	0,1	0,1
Plate Ø	mm	160	160	160	160	160	160	160
Output data	Rs232C	yes	yes	yes	yes	yes	yes	Yes
Dimensions LxWxH	mm	175x220x70	175x220x70	175x220x70	175x220x70	175x220x70	175x220x70	175x220x70
Power supply		Battery DC9V	Battery DC9V	Battery DC9V	Battery DC9V	Battery DC9V	Battery DC9V	Battery DC9V
Item		630.4400.01	630.4600.01	630.4800.01	630.5100.01	630.5120.01	630.5210.01	630.5310.01
Calibration weight	g	200	200	200	500	500	1000	1000

INDUSTRIAL BALANCES


ACS SERIES



Series		ACS 3	ACS 6	ACS 15	ACS 30
Capacity	kg	0 -3	0 -6	0 -15	0 -30
Precision	g	1,0	2,0	5,0	10,0
Plate LxW	mm	245x215	245x215	345x230	345x230
Output data	Rs232C	no	no	no	no
Dimensions LxWxH	mm	250x320x110	250x320x110	350x330x120	350x330x120
Power supply	V/Hz	220/50	220/50	220/50	220/50
Item		630.6003.01	630.6006.01	630.6015.01	630.6030.01

TCS SERIES

Series		TCS 60	TCS 60	TCS 100	TCS 150
Capacity	kg	0 -60	0 -60	0 -100	0 -150
Precision	g	20,0	20,0	30,0	50,0
Plate LxW	mm	300x400	400x500	400x500	400x500
Output data	Rs232C	no	no	no	no
Dimensions LxWxH	mm	300x510x90	400x600x100	400x600x100	400x600x100
Power supply	V/Hz	220/50	220/50	220/50	220/50
Item		630.6060.01	630.6060.02	630.6100.02	630.6150.02



CALIBRATION WEIGHTS

NIKEL PLATED BRASS



Gramms	Item
60 = 50 + 10	630.4000.00
100	630.4001.00
200	630.4002.00
500	630.4005.00
1000	630.4010.00
2000	630.4020.00
5000	630.4050.00

Weighing base equals 205 g.

Composed of:

1 piece 5g. - 1 piece 10 g.

2 pcs. 20g. - 1 piece 50g.

1 piece 100 g

Item	630.4060.00
------	--------------------